

101589220
CA 00/1139

10 NOV 2000 (10.11.00)

PA 312541

REC'D 07 DEC 2000

WIPO

PCT

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

October 10, 2000

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE UNDER 35 USC 111.

APPLICATION NUMBER: 60/156,714

FILING DATE: September 30, 1999

PRIORITY DOCUMENT

SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH RULE 17.1(a) OR (b)



By Authority of the
COMMISSIONER OF PATENTS AND TRADEMARKS

L. Edelen

L. EDELEN
Certifying Officer

Sep-30-99 10:28

From-SIM MCBURNEY

4165951163

T-382 P.03/10 F-557

Please type a plus sign (+) inside this



Docket Number

10024-2JHW

PROVISIONAL APPLICATION FOR PATENT COVER SHEET (Small Entity)

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53 (c).

| INVENTOR(S)/APPLICANT(S) | | | | | |
|--|------------------------|--|---|--------------|------------------------|
| Given Name (first and middle (if any)) | Family Name or Surname | Residence (City and either State or Foreign Country) | | | |
| Michael Gordon L. | THOMPSON HAYWARD | Toronto, Ontario, Canada | | | |
| | | Guelph, Ontario, Canada | | | |
| <input type="checkbox"/> Additional inventors are being named on page 2 attached hereto | | | | | |
| TITLE OF THE INVENTION (280 characters max) | | | | | |
| TRAVERSE SHEAR MODE PIEZOELECTRIC CHEMICAL SENSOR | | | | | |
| CORRESPONDENCE ADDRESS | | | | | |
| Direct all correspondence to: | | | | | |
| <input type="checkbox"/> Customer Number | | <input type="text"/> | | | |
| OR | | <input type="checkbox"/> Place Customer Number Bar Code Label here | | | |
| <input checked="" type="checkbox"/> Firm or Individual Name | | John H. Woodley | | | |
| Address | | Sim & McBurney | | | |
| Address | | 330 University Avenue, 6th Floor | | | |
| City | | Toronto | State | Ontario | ZIP |
| Country | | Canada | Telephone | 416-595-1155 | Fax |
| | | | | M5G 1R7 | 416-595-1163 |
| ENCLOSED APPLICATION PARTS (check all that apply) | | | | | |
| <input checked="" type="checkbox"/> Specification | Number of Pages | 7 | <input type="checkbox"/> Small Entity Statement | | |
| <input type="checkbox"/> Drawing(s) | Number of Sheets | | <input type="checkbox"/> Other (specify) | | |
| METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT (check one) | | | | | |
| <input type="checkbox"/> A check or money order is enclosed to cover the filing fees | | | | | |
| <input checked="" type="checkbox"/> The Commissioner is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number. | | | | | |
| | | | | 182253 | FILING FEE AMOUNT (\$) |
| | | | | | \$75.00 |
| The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government | | | | | |
| <input checked="" type="checkbox"/> No | | | | | |
| <input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are | | | | | |

Respectfully submitted,

SIGNATURE

Date

September 30, 1999

TYPED or PRINTED NAME

John H. Woodley

REGISTRATION NO.
(if appropriate)

27093

TELEPHONE

416-595-1155

USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

SEND TO: Box Provisional Application, Assistant Commissioner for Patents, Washington, DC 20231

TRAVERSE SHEAR MODE

PIEZOELECTRIC CHEMICAL SENSOR

Field of the Invention

5 This invention relates to a process of detecting specific molecules in a liquid (the analyte) with receiving molecules, (the receptors) which are attached to the surface of a thickness shear mode acoustic sensor (TSM). Acoustic energy generated in the sensor is transferred to and from the fluid depending on the surface coupling behaviour. The coupling is altered when the
10 analyte binds to the receptor producing easily measured changes in the electrical characteristics of the sensor.

The invention further relates to the application of the measurement of the coupling effects to the sensing of biomolecules, and other molecules of biological significance such as drugs, in general. For example, the receptor may
15 be a protein, a single oligonucleotide strand, DNA or RNA and the analyte a protein, drug or complementary strands of DNA or RNA. The interaction between the analyte and the sensor bound receptor can be identified through a quantitative TSM response. Other measurement scenarios are possible through the detection of changes in the acoustic coupling between the sensor surface and
20 the surrounding liquid.

Background of the Invention

A TSM sensor is a device which generates mechanical vibrations from an electrical signal and uses these vibrations to detect and/or quantify particular chemical or biochemical substances present in a medium surrounding the
25 sensor (the analyte). Acoustic energy is stored and dissipated both in the device itself, and through interfacial coupling, in a surrounding liquid medium. By coating the sensor with one or more layers of a substance which interacts with the analyte, the energy storage and transfer processes change when the interaction occurs. This changes the acoustic resonance of the sensor, which

60156714-093009

can be observed by measuring the electrical impedance of the sensor. We have published several papers in this field and they are listed as follows:

- 5 1) F. Ferrante, A.L. Kipling and M. Thompson, "Molecular Slip At
The Solid-Liquid Interface Of An Acoustic Wave Sensor", *J.*
Appl. Phys. 76(6):3448-3462, 1994;
- 10 2) G.L. Hayward and M. Thompson, "A Transverse Shear Model Of
A Piezoelectric Chemical Sensor", *Amer. Inst. Physics*
83(40):2194-2201, 1998;
- 15 3) Cavic B.A. et al., "Acoustic Waves And The Real-Time Study Of
Biochemical Macromolecules At The Liquid/Solid Interface",
Faraday Discuss. 107:159-176, 1997;
- 20 4) H. Su and M. Thompson, "Rheological And Interfacial Properties
Of Nucleic Acid Films Studies By Thickness-Shear Mode Sensor
And Network Analysis", *Can. J. Chem.* 74:344-358, 1996.

20 There are several mechanisms whereby a TSM sensor responds to
chemical change on its surface when it is immersed in a liquid. Surface mass
deposition occurs when the analyte binds to the receptor on the sensor surface.
This increases the storage of acoustic energy through the inertia of the added
mass. Acoustic energy may also be stored through the elastic deformation of a
25 coating on the surface. The elasticity of the coating may also change when the
analyte binds to the receptor coating. These energy storage modes determine
the resonant characteristics of the sensor which can easily be measured
electrically. These processes are well known. Examples of piezoelectric
sensors are described, for example in U.S. Patents 5,374,521 and 5,658,732.

30 Viscous loading occurs when acoustic energy is transferred to the liquid.
Some of the acoustic energy is stored by the inertia of the fluid moving with

the sensor surface and can be transferred back to the sensor, but acoustic energy is also dissipated by internal friction within the fluid. The viscous loading effect is also well known, however in the current use of this effect, the transfer of acoustic energy at the surface is considered to be perfect, that is, there is no slip between the sensor surface and the adjacent fluid molecules.

The current practice is based on the well known Butterworth - van Dyke model of a piezoelectric resonator which consists of a resistor, inductor and capacitor in series, all in parallel with another capacitor. The series arm of this network is called the motional arm. Further details of this model and the calculation of their following parameter may be found in the above paper entitled "Rheological and Interfacial Properties of Nucleic Acid Films Studies by Thickness-Shear Mode Sensor and Network Analysis".

Motional Conductance

The motional inductance, L_M , represents the inertial energy stored by the sensor. It depends on the mass of the TSM sensor as well as the mass of material (the analyte) added to the surface. Since liquid coupled to the surface can store and return acoustic energy, L_M is also dependent on the viscosity of the liquid.

Motional Resistance

The motional resistance, R_M , is intrinsically related to the energy dissipated by the sensor.

Accordingly, any imposition of material (or loss of material) that has a viscous property or changes in the viscosity of the liquid will result in a change in the energy dissipation and hence R_M .

Motional Capacitance

The motional capacitance, C_M , represents the elastic energy stored by the sensor. The absorption or chemical binding of the analyte to the coating can have a large effect on the viscoelastic properties of the coating. Depending on

Sep-30-89 10:30 From: SIM MCBURNEY

4165951163

T-382 P.07/10 F-557

the thickness, an added (or removed) layer of material may change the elasticity of the sensor and thus affect C_M . Although most fluids are considered to be viscous, at the high frequencies used in piezoelectric quartz sensors, the liquid may also have elastic properties.

5

Static Capacitance

The static capacitance C_0 represents the dielectric constant of the quartz, but includes that of the medium through the finging field. Charge interactions between the analyte and the sensor coating will affect this value.

10

Summary of the Invention

According to an aspect of the invention, a process for sensing biological or chemical changes in the molecular structure and the mass of the molecules attached to the surface of a transverse shear acoustic device driven by a network analyzer capable of driving the device at a series of sequential predetermined frequencies, measuring electrical parameters at each frequency and fitting a coupling model to these measured parameters, the parameters being F_s , R_M , L_M , C_M and C_0 and parameter α , the process comprises:

- i) exciting the acoustic device at one of the predetermined frequencies for a period of time sufficient to measure the electrical parameters and repeatedly exciting the device at the next sequential frequencies and repeatedly conducting the measurements;
- ii) analyzing the data to determine values for series resonance frequency shift (F_s), motional resistance (R_M), motional inductance (L_M), motional capacitance (C_M), electrostatic capacitance (C_0), mass and boundary slip parameter (α);
- iii) correlating such changes in the parameters with changes in molecular structure and/or mass deposition on the device surface.

Detailed Description of the Preferred Embodiments

This invention is based on the measurement of phenomena based on imperfect acoustic coupling between the sensor surface and the liquid. The

nature of this coupling determines the strength of the viscous loading and elastic effects depending on such parameters as the surface free energy and the molecular conformation of the sensor coating. These molecular parameters are very sensitive to chemical changes at the surface, therefore acoustic coupling

5 provides a novel sensing mechanism.

The impedance measurements are carried out by applying an electrical signal of known frequency and voltage to the sensor and measuring the current through the sensor. Through Ohm's law, this provides the impedance at the known frequency. By performing this measurement over a range of
10 frequencies, a set of data is generated. The above specifically selected parameters of L_M , R_M , C_M and C_O have been found to be the determining parameters for indicating a mass or conformation change at the TSM surface. Hence these parameters are fitted to the data.

While the Butterworth - van Dyke model provides useful information, it
15 is an electrical analogy which presents the information unclearly. An alternate model of the TSM sensor is based on a solution of the equations of motion and electric fields. With this second model as set out in the aforementioned paper entitled "Molecular Slip At The Solid-Liquid Interface Of An Acoustic Wave Sensor" and "A Transverse Shear Model Of A Piezoelectric Chemical
20 Sensor", the deposited mass and the coupling may be determined directly by fitting the electrical impedance data obtained as above. The coupling is represented by a slip parameter, α , which arises from a slip boundary condition used in solving the set of equations. The common approach is to assume perfect coupling and to set $\alpha = 1$. In this invention, α is taken to be a complex number
25 which is determined by fitting the measured impedance data.

The sensing process is understood to be occurring at the interface between the solid device and the liquid medium. Ligands for biological macromolecules include small molecules, ions, proteins, peptides, and strands of both DNA and RNA. The interaction of these entities with the biological
30 molecules attached to the sensor will cause an alteration of the physical properties of the film resulting, in turn, in changes in the measured parameters.

These changes will very clearly result from a combination of some or all of the above response mechanisms particular for each chemical situation. In this regard, the dimensions of the newly bound ligand is an important consideration.

5 The signalling species coated onto the acoustic biosensor are proteins (antibodies, enzymes, hormones, molecular receptors, etc.) and nucleic acids (oligonucleotides, DNA and RNA) attached to the device surface. These molecules exist in a highly hydrated form which can be considered to constitute very viscous lossy gels.

10 The effect of viscous loading is the result of acoustic energy transfer to and from the surrounding medium. This in turn depends on the nature of the contact between the surface and the medium. The contact is controlled by such chemical properties as hydrogen bonding, dispersion interactions and interfacial charge. The process can be viewed as a drag existing between the
15 surface coating and the liquid. α represents the coupling strength but also contains phase shift information. This provides additional information regarding relative mass of liquid molecules compared to those of the sensor surface and when correlated with the selected Butterworth - van Dyke model provide a determination on what is happening at the TSM surface, namely,
20 mass and/or molecular structural shift or change in conformation.

Example

 The human immunodeficiency virus type 1 (HIV-1) is strongly regulated at the transcriptional level by the interaction of an 86-amino acid protein, Tat, with the trans activation responsive element at the 5' -end of the viral
25 messenger RNA transcript (TAR). The TAR-Tat system is an important target for drug discovery research because the binding of the regulatory protein to TAR can be blocked by small molecules.

 In this application we compute the slip parameter α , for the binding of Tat-derived peptides to TAR immobilized on a sensor surface. The TAR RNA
30 is synthesized, with a biotin moiety at the 5' -end, on a DNA synthesizer by standard phosphoramidite chemistry. The acoustic wave sensor is incorporated

into a flow-through configuration and electrically connected to an acoustic network analyzer. A dispersion of 100-500 μ l of the reagent neutravidin is injected into the apparatus and the protein adsorbs to the gold electrode surface of the acoustic wave sensor. A second dispersion of biotinylated TAR- RNA

- 5 (100-500 μ l) is introduced into the system where the formation of the biotin-avidin complex results in attachment of TAR to the sensor surface. Various Tat-derived peptides are then introduced into the flow-through system. In this particular application the following peptides are specified: tat₁₂, tat₂₀, and tat₃₀ where the subscript refers to the number of amino acids in the peptide.
- 10 Dispersions of peptide (100-500 μ l) are injected into the system. On binding of peptide to TAR in real time transient responses in the aforementioned parameters are obtained. The computed α parameter for the various responses, which distinguishes the nature on binding, are as follows:

- 15 Tat₁₂ baseline 1.978 @20.85 degrees
signal 1.964 @ 20.97 degrees

Tat₂₀ baseline 1.985 @21.42 degrees
signal 1.926@ 18.15 degrees

- 20 Tat₃₀ baseline 1.982 @ 22.61 degrees
signal 1.994 @ 23.03 degrees

- 25 Tat₁₂ displays a small decrease in slip magnitude with an increase in phase, whereas tat₂₀ shows large decreases in magnitude and phase. Tat₃₀ depicts smaller increase in magnitude and phase.

Although preferred embodiments of the invention have been described herein in detail, it will be understood by those skilled in the art that variations may be made thereto without departing from the spirit of the invention.

EPM TC 1600

FINAL SEARCH DATE 12/18/03

DELIVER TO GOVT DATE